

**VETERINARY LABORATORY SERVICES UNIT
DEPARTMENT OF VETERINARY LABORATORY DIAGNOSTICS
FACULTY OF VETERINARY MEDICINE, UPM**

SPECIMENS SUBMISSION FOR POST MORTEM LABORATORY					
Type of Specimen	Tests		Specimen Collection	Storage and Transportation	Fee Charges (RM)
Carcass	Complete post mortem	Small sized animal	Carcass	Animals should be kept refrigerated/chilled prior to submission if submitted within 48 hours following death. Animals that would be submitted later than a 48-hours postmortem interval should be frozen	55.20
		Large sized animal	Carcass	Small animals intended for necropsy should arrive packaged either in a carrier, box, or bag. For long transport please include gel ice packs (not wet ice)	237.80